

SWEAR HEALTHCARE

PRODUCT CATALOGUE



choose swear **CHOOSE HEALTH**
Elevate patient care with Swear Healthcare's range of reliable,
high-quality medical disposables



+91-813-082-6865

info@swearhealthcare.com

swearhealthcare.com



SWEAR
HealthCare Pvt. Ltd.

OUR GLOVES

Our wide range of products cater to a variety of medical specialties, providing unparalleled comfort and protection for healthcare professionals and their patients. Discover the peace of mind that comes with wearing Swear Healthcare Gloves.

OUR BRANDS

mYtouch
PRO

mYtouch


medicmile


medicmile⁺Plus



DRUG
LICENSE





STERILE LATEX SURGICAL GLOVES (POWDERED)

Material	: Natural Rubber Latex	Internal Surface	: Powdered
Shape	: Anatomic	Shelf Life	: 5 Years
Cuff	: Beaded	Sterilization	: EO
External Surface	: Textured (Finger-to-palm)	Basic UDI DI	: 89043733LSG4H

SIZE	As Per ISI 13422:1992			As Per Swear Healthcare		
	Min. Length (mm)	Width (mm)	Min. Thickness (mm)	Min. Length (mm)	Width (mm)	Min. Thickness (mm)
6.0	265	76 +/- 6	0.1 for all sizes and areas (Cuff, Palm & Finger)	280	76 +/- 6	0.12 for all sizes and areas (Cuff, Palm & Finger)
6.5	265	83 +/- 6		280	83 +/- 6	
7.0	265	89 +/- 6		280	89 +/- 6	
7.5	265	95 +/- 6		280	95 +/- 6	
8.0	265	102 +/- 6		280	102 +/- 6	

Specifications	As per Standard	As per Swear Healthcare
Tensile Strength (Before Aging)	24 Mpa	24 Mpa
Tensile Strength (After Aging)	18 Mpa	22 Mpa
Ultimate Elongation (%) (Before Aging)	750%	800%
Ultimate Elongation (%) (After Aging)	560%	600%
Stress at 500% Elongation	3.5 Mpa	5.5 Mpa



STERILE LATEX SURGICAL GLOVES (POWDER-FREE)

Material : Natural Rubber Latex

Shape : Anatomic

Cuff : Beaded

External Surface : Textured (Finger-to-palm)

Internal Surface : Polymer Coated

Shelf Life : 5 Years

Sterilization : EO

Basic UDI DI : 89043733LSG4H

SIZE	As Per ISI 13422:1992			As Per Swear Healthcare		
	Min. Length (mm)	Width (mm)	Min. Thickness (mm)	Min. Length (mm)	Width (mm)	Min. Thickness (mm)
6.0	265	76 +/- 6	0.1 for all sizes and areas (Cuff, Palm & Finger)	280	76 +/- 6	0.12 for all sizes and areas (Cuff, Palm & Finger)
6.5	265	83 +/- 6		280	83 +/- 6	
7.0	265	89 +/- 6		280	89 +/- 6	
7.5	265	95 +/- 6		280	95 +/- 6	
8.0	265	102 +/- 6		280	102 +/- 6	

Specifications	As per Standard	As per Swear Healthcare
Tensile Strength (Before Aging)	24 Mpa	24 Mpa
Tensile Strength (After Aging)	18 Mpa	22 Mpa
Ultimate Elongation (%) (Before Aging)	750%	800%
Ultimate Elongation (%) (After Aging)	560%	600%
Stress at 500% Elongation	3.5 Mpa	5.5 Mpa



LATEX EXAMINATION GLOVES (POWDERED)

Material	: Natural Rubber Latex	Internal Surface	: Powdered
Shape	: Ambidextrous	Shelf Life	: 3 Years
Cuff	: Beaded	Sterilization	: EO
External Surface	: Textured	Basic UDI DI	: 89043733LEXGUZ

SIZE	MIN LENGTH (MM)		PALM WIDTH (MM)		WEIGHT PER GLOVE (G)	MINIMUM THICKNESS		
	IS 15354/ ASTM 3578	SWEAR	IS 15354/ ASTM 3578	SWEAR		LOCATION	IS 15354/ ASTM 3578	SWEAR
Small	220	230	80±10	84-85	5.5-6.0	Palm	0.08	0.10
Medium	230	230	95±10	95-96	6.0-6.5			
Large	230	230	111±10	106-107	6.5-7.0	Finger	0.08	0.11
X-Large	230	230	120±10	112-113	7.0-7.5			

Properties	Before Ageing		After Ageing (70 ± 2° C for 166 ± 2 h)	
	IS 15354/ ASTM 3578	SWEAR	IS 15354/ ASTM 3578	SWEAR
Tensile Strength (mpa)	18 min	20 min	14 min	18 min
Ultimate Elongation (%)	650 min	800 min	500 min	600 min
Stress at 500% Elongation (Mpa)	5.5 max	5.5 max	N/A	N/A

ANTI-MICROBIAL SURGICAL GLOVES (POWDERED & POWDER FREE)

Our Anti-Microbial Examination Gloves (Powdered and Powder-Free) are perfect for medical procedures. These ambidextrous gloves feature a comfortable rolled rim and are safe for handling contaminated materials. Providing 99.9% protection against microbial infections, they ensure a smooth, low-allergen experience. Our gloves conform to ASTM D7909, ASTM D3577, EN455, and ASTM D6124 standards for maximum safety and reliability.

Details:



6, 6.5, 7, 7.5 AND 8



OUR GLOVE'S ANATOMY

285 MM EXTENDED LENGTH
offer greater coverage than
standard options.

MICRO-ROUGHENED TEXTURE

from finger tips to entire palm for
improved dexterity.



MEDICAL-GRADE PAPER POUCHES
ensure superior packaging with
optimal glove sterility and protection.






UROMETER

URINE COLLECTION BAG WITH WITH
MEASURED VOLUME METER

NON TOXIC | NON PYROGENIC | STERILE | SINGLE USE


- **Precision in Monitoring:** Designed for accurate measurement of both small and large urine output with a clear, graduated **250 ml** reservoir attached to a **2000 ml** overflow bag.
- **Enhanced Safety:** The closed drainage system and "T"-shaped bottom outlet reduce the risk of contamination and catheter-associated urinary tract infections (CAUTIs).
- **User-Friendly Design:** Equipped with **120 cm kink-resistant tubing** for efficient drainage and a foldable hanger with a tube holder for convenient bedside attachment.



URINE COLLECTION BAG

WITH TOP OUTLET & NON RETURN VALVE
NON-TOXIC | LATEX FREE | SINGLE-USE

- **Enhanced Transparency:** Made from a **300-micron** clear and transparent PVC sheet, ensuring easy visibility to detect any colour changes or contamination.
- **Hygienic Drainage Valve:** Features a **top conical drainage** valve that allows for efficient bag emptying while minimizing direct contact with urine, ensuring better hygiene and safety.
- **Kink-Free Tubing:** Equipped with kink-free tubing of **100 cm length** that remains durable and effective for ensuring uninterrupted flow



100% automated manufacturing unit

to deliver precise, efficient, and high-quality sustainable gloves that meet the needs of healthcare professionals and patients.

Why choose us?

Advanced Technology: 100% automatic production lines for precision and efficiency.

Global Reach: Exporting to 28+ countries all around the World.

Sustainable Practices: Committed to producing the highest quality gloves using sustainable methods.

CONTACT



335 & 336 (3rd floor) Paras Trade Center,
Gwal Pahadi, Gurgaon - 122003, India



+91-813-082-6865



swearhealthcare.com



RIICO Growth Center (Extension, F-374)
Dholpur, Sanda, Rajasthan - 328001, India




customercare@swearhealthcare.com
info@swearhealthcare.com



Currently exporting to 28+ countries

With a strong focus on innovation and quality, Swear Healthcare has positioned itself as a leader in the global market for latex gloves.



CORPORATE OFFICE

Unit No. 335 & 336 (3rd Floor)
Paras Trade Center, Gwal Pahadi
Gurugram, Haryana, India - 122003

MANUFACTURING UNIT

F-374, RIICO Growth Center
Extension, Dholpur, Rajasthan,
India - 328001